

# SL-D701 Digital Slit Lamp

High-resolution imaging and video for accurate documentation and effective patient education

- ✓ Modern slit lamp ergonomics + digital imaging with the DC-4 camera
- ✓ Designed for advanced ocular surface and dry eye evaluation



LIMITED-TIME OFFER


# \$5,000

instant credit when you upgrade to Digital Slit lamp with DC-4 camera\*

DC-4 Digital  
Camera Attachment

## ENHANCED FLUORESCEIN EVALUATION

Optional yellow filter highlights fluorescein for clearer evaluation of ocular surface staining and contact lens fit.



Yellow + Cobalt Blue Filter<sup>1</sup> Cobalt Blue Filter<sup>2</sup>

1. Image courtesy of Paymaan Asnaashari, OD, FAAO, FSLs. Arden Park Optometry
2. Image courtesy of Toru Noda, M.D. Department of Ophthalmology, Tokyo Medical Center

## MEIBOMIAN GLAND IMAGING

Infrared illumination enables detailed visualization and documentation of meibomian gland structure and health.



Meibomian gland imaging captured with SL-D701 slit lamp + DC-4 Digital Camera + BG-5 background illuminator.<sup>2</sup>

2. Image courtesy of Reiko Arita, MD, PhD, Itoh Clinic, Department of Ophthalmology, University of Tokyo School of Medicine.

\*Offer not valid for accounts participating in buying group or contracted pricing programs.